

ENG **ULTRASONIC NEBULIZER LD.** Instruction Manual.
POL **INHALATOR ULTRADZIĘKIOWY LD.** Instrukcja obsługi.
HUN **ULTRAHANGOS LD INHALÁTOR.** Használati utasítás.
ROU **INHALATOR CU ULTRASUNETE LD.** Manual de instrucțiuni.
BGR **УЛТРАЗВУКОВ ИХАЛАТОР LD.** Ръководство за експлоатация.

fig.1 rys. 1 **1. ábra** **имaginea 1** **фиг.1**
PARTS AND COMPONENTS **PARTI DE BAZĂ A DISPOZITIVULUI** **ОСНОВНИ ЧАСТИ НА УРЕДА**
ПОДСТАВОВЕ ЧІСЛІ КОМПОНЕНТЫ





fig.2, 2b, 2c rys. 2, 2b, 2c **2a, 2b, 2c ábra** **имaginea 2, 2b, 2c** **фиг.2, 2b, 2c**
OPERATING PROCEDURES **SPOSÓB UŻYTKOWANIA** **ПОСЛЕДОВАТЕЛЬНОСТЬ НА УПОРЯДОЧЕНИЕ**





fig.3 rys. 3 **3. ábra** **имaginea 3** **фиг.3**

THE PARTICLES OF DIFFERENT SIZES ARE DISTRIBUTED AEROSOLS AS FOLLOWS*:
AEROZOŁO ROZNEJ WIELKOŚCI CZĄSTECZEK ROZPROWADZANY JEST W NASTĘPUJĄCY SPOSOB:
A KÜLNÖBÖZŐ MÉRÜT RÉSZCSEKK AEROSZOL ELOSZLÁSA A KÖVETKEZŐ:
PARTICULES DE DIVERSE DIMENSIUNI SE DISTRIBUIE ÎN AEROSOL ÎN FELUL URMĂTOR*:
ЧАСТИЦЕ СРАЗЛИЧЕН РАЗМЕР СЕ РАЗПРЕДЕЛЯТ В АЕРОЗОЛА ПО СЛЕДНИЯ НАЧИН*:



ENG
PARTS AND COMPONENTS (fig.1)

1. Nebulizer chamber
2. Nebulizer chamber
3. The upper part of the nebulizer chamber
4. Nebulizer chamber cap
5. The middle part of the nebulizer chamber
6. Gasket
7. Solution container
8. The lower part of the nebulizer chamber
9. Main unit
10. Piezoelectric element
11. Retainer indicator
12. Air fan
13. Fuse 250V/1,5A
14. Fuse 250V/0,5A
15. Power plug
16. Adult mask
17. Nebulizer solution container
18. Nebulizer mouthpiece
19. Inhalation tube
20. Joint sleeve
21. Bag

GENERAL INFORMATION
 Ultrasonic nebulizer (hereinafter referred to as nebulizer, LD-152U, device) is intended for the treatment and prevention of respiratory diseases with aerosols of AQUEOUS drug solutions in hospitals and at home. This manual is intended to assist the user in the safe and efficient operation of Ultrasonic Nebulizer LD-152U. The instrument should be used in accordance with the rules stated in this manual and should not be used for the purposes other than those described herein. It is important to read and understand the whole manual. Functionally, the device consists of an aerosol generation chamber (hereinafter: nebulizer chamber), a piezoelectric element, and a control circuit. A fan is located in the body for cooling the electronic components and forcing air into the nebulizer chamber. On the body of the device, there are on/off controls, regulators for the intensity of aerosol formation, and an airflow regulator.

CONTRAINDICATIONS
 Major neonatal, systemic, blood diseases, severe general exhaustion, stage III hypertension, pronounced arteriosclerosis or the ventricular fibrillation of the cardiovascular system in the stage of decompensation, kwashiorrhea and sun stroke, also the presence of the ventricular fibrillation of the cardiovascular system in the stage of decompensation, kwashiorrhea and sun stroke problems. Allergy to the nebulizer chamber, the nebulizer mouthpiece, the inhalation tube, the nebulizer solution container, the piezoelectric element and the control circuit. A fan is located in the body for cooling the electronic components and forcing air into the nebulizer chamber. On the body of the device, there are on/off controls, regulators for the intensity of aerosol formation, and an airflow regulator.

PRECAUTIONS
Attention! It is not allowed to use nebulizer solutions containing ether, oils, or suspended particles (suspensions), including decoctions and herbal infusions. All types of standard nebulizer solutions in liquid form for nebulizer therapy produced by pharmaceutical companies are recommended for use.

Attention! This device is not intended for inhalation with solutions of hormonal drugs and antibiotics, because the use of the piezoelectric element leads to their destruction.

After opening the package, check the integrity of the device, paying particular attention to the presence of damage to plastic parts, which can make accessible parts of the under-voltage device, as well as damage to the power cable. Switch on the device only after making sure that there is no such damage. Perform this check before each use;

2. Before connecting the device, make sure that the electrical characteristics of the device correspond to those of the mains supply.
 3. When you are not using the device, disconnect it from the mains by pulling out the plug from the socket.
 4. Observe the safety regulations for electrical appliances, in particular:
 - Never immerse the device in water;
 - Install the device vertically on a flat, stable surface;
 - Position the device so as not to block the ventilation holes on the bottom of the device;
 - Do not use this device in explosive environments, e.g., near flammable anesthetics or inside an oxygen chamber;
 - Do not touch the power plug with wet hands;
 - Do not touch the power cord while in water or steam;
 - Do not touch the device if it has fallen into a salt bath. Do not use when bathing;
 - The device and especially its accessories should be kept out of the reach of children, as they contain small parts that can be swallowed;
 - Do not pull on the power cord to remove the plug from the socket;
 - Do not touch the power plug with wet hands;
 - Do not store the device indoors, away from heat sources. After use, it is recommended to store the device away from dust and direct sunlight;

5. If the device does not work, see the section "TROUBLESHOOTING ALGORITHM AND REMEDIES".

6. The device is not intended for outdoor use.
 7. None of the electrical or mechanical parts are intended to be repaired by end users. Do not open the device yourself. For repairs, contact exclusively the authorized service centers listed on the website www.littledoctor.sg
Technical data
Model LD-152U
Type ultrasonic
Power consumption, no more 40 VA
Aerosol production performance 2 ml/min.*
Mean median diameter of aerosol particles (MMAD) 3 µm*
Ultrasound frequency 1.7 MHz ± 10 %
Water temperature < 60°C
Volume of the nebulizer solution container 12 mL
Residual volume of nebulizer solution, no more 1 mL
Noise level, no more 55 dB*
Time of continuous work up to 4 hours
Cooldown time of the device at least 40 min
Power supply: 230V 50Hz
Protection against electric shock Type B
Device operational conditions:
 - Ambient temperature: 10°C to 40°C
 - Humidity: 15% to 85% RH
Storage and transportation conditions for the device:
 - Ambient temperature: -20°C to 40°C
 - Humidity: 15% to 85% RH
Weight of the device (without packaging) 1,515 kg
Overall dimensions of the electronic unit 250x104x160 mm
Indicated on the body of the device in the serial number as "AYMM15XXXXX", where YY is the year, and MM is the month of manufacture.

* - data obtained by Little Doctor International (S) Pte. Ltd.
 This device is certified according to international standard ISO 13485. Device comply with the requirements of European Directive MDD 93/42/EEC.

SYMBOL EXPLANATION
 2. Conditions for storage, transportation and use
 4. Representative in the European Union

The revision date of this Operation Manual is shown on the last page as EXXX/YYMM/MM, where YY is the year, MM is the revision month, and the MM is the revision number.
TROUBLESHOOTING ALGORITHM AND REMEDIES
 In the absence of aerosol formation, troubleshooting is recommended according to the following order:

OPERATING PROCEDURES
Important! Before using the device for the first time, it must be thoroughly cleaned, as described in paragraph 1 of the "Maintenance, storage, repair, and disposal" section

1. Remove the nebulizer chamber from the main unit of the device, holding it by the tabs on the nebulizer chamber cap (Figure 2a).
 2. Fill the container with water exactly up to the "WATER LEVEL" mark on the level sensor. If the water level is less than the permissible level, the sensor will not allow the piezoelectric element to be turned on, while the red indicator will be on (Figure 2b).
 3. Place the nebulizer chamber vertically on a counter-clockwise. So that the marker « » points to « » (Figure 2b). Pull the cap of the nebulizer chamber up and remove it.

Attention! Do not turn on the device with an empty drug container. Do not add more than 12 mL of nebulizer solution to the container.

Attention! Do not turn on the device with an empty drug container so that the marker « » points to « » (Figure 2c).

5. Close the nebulizer chamber cap by turning it clockwise so that the marker « » points to « » (Figure 2c).
 6. Place the nebulizer chamber back into the main unit of the device.
 7. Connect the tube with a mask or mouthpiece to the nebulizer chamber. If necessary, you can lengthen the tube by stretching it or connecting the two tubes with a sleeve.

It should be born in mind that increasing the length of the inhalation tubes leads to an increase in the loss of the nebulizer solution.

Important! Each patient is advised to use an individual mouthpiece and/or mask.

8. Set the controls to their original positions:
 - regulator « » to position « »;
 - airflow regulator « » to position « »;
 - power plug « » to position « »

Important! For effective inhalation, you need to calm down, relax, and sit up straight.

If you inhale in bed, place a pillow under your back so that your back is as straight as possible. An awkward posture and uneven breathing can block some parts of the lungs and the airway.

10. Start inhalation, turn the regulator « » clockwise to the position "1", the green indicator « » will light up, and the device will begin to produce an aerosol.

12. In the initial period of inhalation, adjust the optimal level of aerosol formation using the regulator « ».

The intensity of aerosol formation, and consequently, the nebulizer solution flow rate depends on the viscosity of the solution used. For example, when using 0.9% sodium chloride as a nebulizer solution, the dependence of the nebulizer solution flow rate on the regulator position « » is shown in the table.

Authorized Representative in the EU: Little Doctor Europe Sp. z o.o., 576 Zawila Street Krakow 30-390 Poland.

For more information please visit www.littledoctor.sg

POL
PODSTAWOWE CZEŚCI I KOMPONENTY (rys. 1)

Regulator position	The approximate nebulizer solution flow time (volume 10 mL, the remainder is not more than 1 mL, the airflow regulator position is «»), min
1	15
2	8
3	6

NOTE 1: If you are using a nebulizer for children, too much aerosol flow can make it difficult for the baby to breathe. In this case, it is recommended to reduce the aerosol flow with the airflow regulator, which will increase the inhalation time.

NOTE 2: During inhalation, condensation may form under the body of the device.

13. After using:
 a) turn off the device by setting the regulator « » to the position "0";
 b) remove the power plug from the socket;
 c) disconnect the nebulizer chamber from the main unit, remove the upper cap of the nebulizer chamber and drain the remaining nebulizer solution;
 d) drain the water from the container;
 e) wipe the device with a soft cloth and let it dry;
 f) move the device to a dry place out of reach of children.

REPLACING THE CONTAINER OF SOLUTION

Inspect the container of the solution before each use. In case of violation of integrity, deformation or other damage, replace the container with a new one.

Proceed as follows to replace the container:

1. Remove the nebulizer chamber from the main unit (pos. 1 and 2 of the "PARTS AND COMPONENTS" section).

2. Unscrew the lower part of the nebulizer chamber (pos. 1).

3. First remove the gasket (pos. 1.4) and then the container (pos. 1.5) from the bottom of the nebulizer chamber.

4. Place the new container in the bottom of the nebulizer chamber and install the gasket on top. Make sure the container and gasket fit tight against the walls.

5. Screw the lower part of the nebulizer chamber to the middle part of the chamber (pos. 1.3).

6. Install the nebulizer chamber on the main unit.

MAINTENANCE, STORAGE, REPAIR, AND DISPOSAL

1. Clean the device and all accessories regularly. Disinfect device accessories after each use. All device accessories are recommended to be washed by water, stirring, rinsing, then autoclaving or sterilization at 120°C for 10 minutes, or autoclaved at 150°C. After treatment, wipe all parts of the device with a soft cloth.

2. At home, using 3-6% alcohol solution. The disinfection time is 30 min. The initial temperature of the solution is 50 °C. The shelf-life of the working solution is 1 day / using a 3-6% solution of table vinegar. The solution is prepared from 1 part of 3-6% table vinegar and 3 parts of distilled water. Make enough solution to immerse the parts.

The disinfection time is 30 min. Do not reuse working solutions! After treatment, rinse all the device accessories under an abundant stream of hot water. Air dry the accessories or blot with a soft, clean, lint-free cloth.

3. Do not expose airways to strong sunlight or shock.

4. Do not store or use the device in the vicinity of heaters and open flames.

5. Protect the device from contamination.

6. Do not allow the device to come into contact with aggressive solutions.

7. If necessary, carry out repairs only with specialized contractors.

8. The life cycle of the drug container is 6 months from the beginning of the operation. The manufacturer has not established special conditions of disposal.

PRECAUTIONS

Attention! It is not allowed to use nebulizer solutions containing ether, oils, or suspended particles (suspensions), including decoctions and herbal infusions. All types of standard nebulizer solutions in liquid form for nebulizer therapy produced by pharmaceutical companies are recommended for use.

Attention! This device is not intended for inhalation with solutions of hormonal drugs and antibiotics, because the use of the piezoelectric element leads to their destruction.

After opening the package, check the integrity of the device, paying particular attention to the presence of damage to plastic parts, which can make accessible parts of the under-voltage device, as well as damage to the power cable. Switch on the device only after making sure that there is no such damage. Perform this check before each use;

TECHNICAL SPECIFICATIONS

Model	LD-152U
Type	ultrasonic
Power consumption, no more	40 VA
Aerosol production performance	2 ml/min.*
Mean median diameter of aerosol particles (MMAD)	3 µm*
Ultrasound frequency	1.7 MHz ± 10 %
Water temperature	< 60°C
Volume of the nebulizer solution container	12 mL
Residual volume of nebulizer solution, no more	1 mL
Noise level, no more	55 dB*

