# LD-300, LD-301, LD-302

**ENG** LD DIGITAL MEDICAL THERMOMETER. INSTRUCTION MANUAL

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
## ENG

Using this LD thermometer you can promptly, accurately and safely measure your body temperature.

The thermometer is intended for measuring body temperature **ONLY** by using the methods specified in this Instruction manual.

Temperature measurement range, °C	32.0 - 42.0	Automatic cutoff	In 10 minutes after measurement is complete
Temperature display range, °C	32.0 - 42.0	Overall dimensions: Dimension (w/o a case), mm	129x18x10 (LD-300), 129x18x11 (LD-301), 129x19x11 (LD-302)
Measuring method	Upper limit rating method	Weight (w/o a case), g	10 (LD-300), 8 (LD-301), 11 (LD-302)
Limits of absolute temperature measurement error, °C	±0.1	Complete set	Thermometer with a battery, plastic case, service manual, package
Probe tip	Thermistor	Advisable calibration interval	Maximum 2 years
Indication	Liquid crystal display with indication to one-tenth of one degree	Year of production	See LOT code (the first 2 numerals) on the back side of the thermometer body.
Units	Degrees Celsius (°C)	Operation conditions: ambient air temperature, °C	From +10 to +35
Beep	Measurement is complete	relative humidity, %Rh	15 - 85
Memory	Last measured value	Storage conditions: ambient air temperature, °C	From -25 to +55
Supply voltage, V	1,5	relative humidity, %Rh	15 - 85
Power source	One LR41 battery		
Maximum power consumption, MW	0,15		
Battery life time	About 100 hours of measurement		

This unit manufacturing is certified according to international standard ISO 13485. Device comply with the requirements of European Directive MDD 93/42/EEC.

### SAFETY MEASURES

- For a waterproof part of the LD-300 thermometer, see Figure (2). Prevent ingress of water in other parts of the device. LD-301 and LD-302 thermometers have a waterproof body that may be submerged in water down to the depth of 16 cm for max 30 minutes.
- This thermometer is designed to measure underarm, oral or rectal temperatures.
- Keep children away of small parts. These parts may be hazardous to health. Use them under control only!

### MEASUREMENT PROCEDURE. MEMORY

- Press a button (4) to switch on the device.
- On pressing the button all segments will be indicated on a display (3), as demonstrated in the figure on the right. Now a thermometer is in active status.
- Release the button. The last temperature reading will be displayed. Then symbols “L °C” will be displayed with flashing symbol “°C”.
- Take the thermometer under the arm by placing the tip into the hollow of the armpit (or in other place – see Section “Temperature measurement methods”.
- If body temperature exceeds 32.0 °C, a temperature value will be displayed in place of symbols “L °C” with a marker “°C” flashing for all a measurement period.
- If temperature values are not changed for 16 seconds, the marker “°C” will stop flashing and a beep signal will sound. Time of measurement depends on a selected measurement method (see Section “Temperature measurement methods”). If the procedure meets the requirements specified in Section “Temperature measurement methods”, a beep signal will indicate to completion of measurement. If you doubt the accuracy of measurement, extend measurement time by 2-3 min after a beep signal sounds. If the measured temperature is below 32 °C, symbols “L °C” will be steadily displayed. If the measured temperature is above 32 °C, symbols “H °C” will be displayed.
- Temperature readings will be displayed till the thermometer is switched off.
- Press the button (4) to switch off the thermometer.
- If you forget to turn off the thermometer, it will turn off automatically in about 10 minutes.

### TEMPERATURE MEASUREMENT METHODS

#### ORAL MEASUREMENT

Insert the thermometer probe into a hygienic packet (option). Before taking temperature keep the mouth closed for 5 minutes. Then place the probe under the tongue with the tip rested in. Close the mouth. Oral temperature shall be taken for 1 minute.

#### RECTAL MEASUREMENT

Gently slide the thermometer tip into the rectum. Make sure that the tip is inserted into the rectum by no more than 2.5 cm. Rectal temperature shall be taken for 1 minute.

#### ARMPIT MEASUREMENT

Before taking temperature readings wipe the hollow of the armpit by a dry towel and keep the arm in down position for 5 minutes. Place the tip well into the underarm while pressing the arm tight. Underarm temperature shall be taken for 3-4 minutes.

NOTE! Rectal temperature should be 1 degree higher than the oral one and underarm temperature – 1 degree less than the oral one. For any questions consult your physician.

### BATTERY REPLACEMENT

If a battery is discharged (less than 1.4V), the symbol «⊞», will be displayed thereby indicating for the need of battery replacement.

- Open the battery compartment (5) sliding the cover opposite to the thermometer probe.
- Remove a discharged battery. Do not use any sharp metal objects.
- Insert a new battery making sure that the positive (+) pole is facing up and the negative (-) pole - down.

### WARRANTY

- The warranty for the device is 24 months after the purchase date.
- The warranty liabilities are contained in the warranty card given at the sale of this Device to a purchaser.
- The addresses of organizations providing warranty service can be found on the website www.littledoctor.sg.

☒ QUALITY CLAIMS ARE RECEIVED AT THE FOLLOWING ADDRESS:

Little Doctor Europe Sp. z o.o. 57G Zawila Street, 30-390, Krakow, Poland. Service phone: +48 12 2684748, 2684749

### SYMBOLS

☐ CE marking in conformity with EC directive 93/42/EEC	☔ Keep dry	☐ LOT	LOT number	⚠ Important: Read the manual
☐ EC Representative	⌚ When utilizing the waste, refer to current rules applicable in your region	☐ SN	Serial number	IP22 IP protection class
☠ Type BF	☞ Storage/operation temperature and humidity	☐	Mark of conformity of the Ukraine	
☐ Manufacturer				

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